



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

FZD8 Antibody [#G261]

Catalog Number: G261

Concentration: 1mg/ml

Swiss-Prot No.: Q9H461

Other Names: Frizzled-8 [Precursor]; Fz-8; hFz8; FZD8

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FZD8.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FZD8 antibody detects endogenous levels of total FZD8 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Saitoh T., Int. J. Oncol. 18:991-996(2001).

Deloukas P., Nature 429:375-381(2004).

