



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ADA2L Antibody [#C10181]

Catalog Number: C10181

Concentration: 1mg/ml

Swiss-Prot No.: O75478

Other Names: Transcriptional adapter 2-like, ADA2-like protein, TADA2L, KL04P

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ADA2L.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ADA2L antibody detects endogenous levels of total ADA2L protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Candau R., Mol. Cell. Biol. 16:593-602(1996).

Ogryzko V.V., Cell 94:35-44(1998).

Henriksson A., Mol. Cell. Biol. 17:3065-3073(1997).

