



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

FZD9 Antibody [#G262]

Catalog Number: G262

Concentration: 1mg/ml

Swiss-Prot No.: O00144

Other Names: Frizzled-9; Fz-9; hFz9; FzE6; CD349 antigen; FZD9; FZD3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human FZD9.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FZD9 antibody detects endogenous levels of total FZD9 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%);
Rat (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Wang Y.-K., Hum. Mol. Genet. 6:465-472(1997).

Hillier L.W., Nature 424:157-164(2003).

Tanaka S., Proc. Natl. Acad. Sci. U.S.A. 95:10164-10169(1998).

