



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TAF5L Antibody [#C10185]

Catalog Number: C10185

Concentration: 1mg/ml

Swiss-Prot No.: O75529

Other Names: TAF5-like RNA polymerase II
p300/CBP-associated factor-associated factor 65 kDa
subunit 5L, PCAF-associated factor 65 beta, PAF65-beta,
TAF5L, PAF65B

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from internal of human

TAF5L.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: TAF5L antibody detects endogenous levels
of total TAF5L protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Ogryzko V.V., Cell 94:35-44(1998).

Stanchi F., Yeast 18:69-80(2001).

Struhl K., Cell 94:1-4(1998).

