



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## FOXH1 Antibody [#C10188]

**Catalog Number:** C10188

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O75593

**Other Names:** Forkhead box protein H1, Forkhead activin signal transducer 1, Fast-1, hFAST-1, Forkhead activin signal transducer 2, Fast-2, FOXH1, FAST1, FAST2

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

FOXH1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** FOXH1 antibody detects endogenous levels of total FOXH1 protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%);  
Mouse (Identities = 100%, Positives = 100%)

**Applications:** IF: 1:100~1:500 ELISA: 1:5000

**References:**

Zhou S., Mol. Cell 2:121-127(1998).

