



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AKAP3 Antibody [#C10200]

Catalog Number: C10200

Concentration: 1mg/ml

Swiss-Prot No.: O75969

Other Names: AKAP110, FSP95, PRKA3, SOB1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human AKAP3.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AKAP3 antibody detects endogenous levels of total AKAP3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:10000

References:

Lefevre A., Biochem. Biophys. Res. Commun. 259:60-66(1999).

Mandal A., Biol. Reprod. 61:1184-1197(1999).

Vijayaraghavan S., Mol. Endocrinol. 13:705-717(1999).

