



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## FIR Antibody [#C10207]

**Catalog Number:** C10207

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O94887

**Other Names:** FERM, RhoGEF and pleckstrin domain protein 2, FGD1-related Cdc42-GEF, FIR, FRG, PLEKHC3

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from Internal of human FIR.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** FIR antibody detects endogenous levels of total FIR protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Nagase T., DNA Res. 5:31-39(1998).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Jandrig B., Oncogene 23:9295-9302(2004).

