



Web: www.anbio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

JHD3B Antibody [#C10209]

Catalog Number: C10209

Concentration: 1mg/ml

Swiss-Prot No.: O94953

Other Names: EC 1.14.11.-, JMJ2B, Jumonji domain containing protein 2B

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human JHD3B.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: JHD3B antibody detects endogenous levels of total JHD3B protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:1000

References:

Nagase T., DNA Res. 5:355-364(1998).

Nakajima D., DNA Res. 9:99-106(2002).

Grimwood J., Nature 428:529-535(2004).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using JHD3B antibody.