



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ALX3 Antibody [#C10213]

Catalog Number: C10213

Concentration: 1mg/ml

Swiss-Prot No.: O95076

Other Names: Homeobox protein aristaless-like 3,
Proline-rich transcription factor ALX3, ALX3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ALX3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: ALX3 antibody detects endogenous levels of total ALX3 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:5000

References:

Wimmer K., Genes Chromosomes Cancer
33:285-294(2002).



Western blot analysis of extracts from 293 cells, using ALX3 antibody.



Immunofluorescence analysis of HeLa cells, using ALX3 antibody.