



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## TNFSF15 Antibody [#C10215]

**Catalog Number:** C10215

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O95150

**Other Names:** Tumor necrosis factor ligand superfamily member 15, Vascular endothelial cell growth inhibitor, TNF ligand-related molecule 1, TNFSF15, TL1, VEGI

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from C-terminal of human

TNFSF15.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


**Specificity:** TNFSF15 antibody detects endogenous levels of total TNFSF15 protein.

**Reactivity:** Human


**Applications:** WB: 1:500~1:3000 IF: 1:100~1:500  
ELISA: 1:20000

**References:**

Zhai Y., FASEB J. 13:181-189(1999).



Western blot analysis of extracts from COLO205 cells, using TNFSF15 antibody.



Immunofluorescence analysis of HeLa cells, using TNFSF15 antibody.