



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GALR1 Antibody [#G267]

Catalog Number: G267

Concentration: 1mg/ml

Swiss-Prot No.: P47211

Other Names: Galanin receptor type 1; GAL1-R; GALR1;
GALR1; GALNR, GALNR1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GALR1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GALR1 antibody detects endogenous levels of total GALR1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 81%, Positives = 87%)

Applications: IF: 1:100~1:500 ELISA: 1:5000

References:

Habert-Ortoli E., Proc. Natl. Acad. Sci. U.S.A. 91:9780-9783(1994).

Lorimer D.D., Biochem. Biophys. Res. Commun. 241:558-564(1997).

Jacoby A.S., Genomics 45:496-508(1997).

