



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** +86 519 8805 0026

## GNRHR Antibody [#G276]

**Catalog Number:** G276

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P30968

**Other Names:** Gonadotropin-releasing hormone receptor; GnRH receptor; GnRH-R; GNRHR; GRHR

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human GNRHR.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** GNRHR antibody detects endogenous levels of total GNRHR protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%);  
Mouse (Identities = 85%, Positives = 100%)

**Applications:** IF: 1:100~1:500 ELISA: 1:1000

**References:**

Kakar S.S., Biochem. Biophys. Res. Commun. 189:289-295(1992).

Chi L., Mol. Cell. Endocrinol. 91:R1-R6(1993) Kakar S.S., Mol. Cell. Endocrinol. 106:145-149(1994).

