



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR115 Antibody [#G285]

Catalog Number: G285

Concentration: 1mg/ml

Swiss-Prot No.: Q8IZF3

Other Names: FLJ38076; G-protein coupled receptor 115;
G-protein coupled receptor PGR18; GP115; PGR18

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GPR115.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR115 antibody detects endogenous levels of total GPR115 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Mungall A.J., Nature 425:805-811(2003).

Vassilatis D.K., Proc. Natl. Acad. Sci. U.S.A. 100:4903-4908(2003).

Fredriksson R., FEBS Lett. 531:407-414(2002).

