



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MED26 Antibody [#C10228]

Catalog Number: C10228

Concentration: 1mg/ml

Swiss-Prot No.: O95402

Other Names: ARC70, CRSP complex 7, CRSP70, MED26, activator-recruited cofactor 70 kDa component, cofactor required for Sp1 transcriptional activation 7

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human

MED26.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MED26 antibody detects endogenous levels of total MED26 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Ryu S., Nature 397:446-450(1999).

Naeae A.M., Nature 398:828-832(1999).

Sato S., Mol. Cell 14:685-691(2004).

