



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

SIX6 Antibody [#C10233]

Catalog Number: C10233

Concentration: 1mg/ml

Swiss-Prot No.: O95475

Other Names: Homeobox protein SIX6, Sine oculis homeobox homolog 6, Optic homeobox 2, Homeodomain protein OPTX2, SIX6, OPTX2, SIX9

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human SIX6.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SIX6 antibody detects endogenous levels of total SIX6 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:5000

References:


Lopez-Rios J., Mech. Dev. 83:155-159(1999).

Gallardo M.E., Genomics 61:82-91(1999).

Gallardo M.E., Am. J. Med. Genet. A 129:92-94(2004).



Western blot analysis of extracts from HUVEC cells and COLO205 cells, using SIX6 antibody.



Immunofluorescence analysis of HepG2 cells, using SIX6 antibody.