



**Web:** [www.anbio.com](http://www.anbio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## CBX6 Antibody [#C10234]

**Catalog Number:** C10234

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O95503

**Other Names:** Chromobox protein homolog 6, CBX6

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human CBX6.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** CBX6 antibody detects endogenous levels of total CBX6 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Dunham I., Nature 402:489-495(1999).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

