



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

GAS41 Antibody [#C10237]

Catalog Number: C10237

Concentration: 1mg/ml

Swiss-Prot No.: O95619

Other Names: YEATS domain-containing protein 4, Glioma-amplified sequence 41, Gas41, NuMA-binding protein 1, NuBI-1, NuBI1, YEATS4, GAS41

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

GAS41.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GAS41 antibody detects endogenous levels of total GAS41 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Fischer U., Hum. Mol. Genet. 6:1817-1822(1997).

Harborth J., J. Biol. Chem. 275:31979-31985(2000).

Ota T., Nat. Genet. 36:40-45(2004).

