



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

ZNF682 Antibody [#C10241]

Catalog Number: C10241

Concentration: 1mg/ml

Swiss-Prot No.: O95780

Other Names: Zinc finger protein 682, ZNF682

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ZNF682.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ZNF682 antibody detects endogenous levels of total ZNF682 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Ota T., Nat. Genet. 36:40-45(2004).

Grimwood J., Nature 428:529-535(2004).

