



Web: [www.anbio.com](http://www.anbio.com)

[www.dijibio.com](http://www.dijibio.com)

Order: [order@dijibio.com](mailto:order@dijibio.com)

Tel: (86) 519 8805 0026

## MBD3 Antibody [#C10250]

**Catalog Number:** C10250

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O95983

**Other Names:** Methyl-CpG-binding domain protein 3, Methyl-CpG-binding protein MBD3, MBD3

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from C-terminal of human MBD3.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** MBD3 antibody detects endogenous levels of total MBD3 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Hendrich B., Mol. Cell. Biol. 18:6538-6547(1998).

Grimwood J., Nature 428:529-535(2004).

Tatematsu K., Genes Cells 5:677-688(2000).

