



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## TOP3B Antibody [#C10251]

**Catalog Number:** C10251

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O95985

**Other Names:** DNA topoisomerase 3-beta-1, EC 5.99.1.2,  
DNA topoisomerase III beta-1, TOP3B, TOP3B1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human TOP3B.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** TOP3B antibody detects endogenous levels of total TOP3B protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Ng S.-W., Nucleic Acids Res. 27:993-1000(1999).

Kawasaki K., Genome Res. 7:250-261(1997).

Collins J.E., Genome Biol. 5:RESEARCH84.1-RESEARCH84.11(2004).

