



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

HAND1 Antibody [#C10253]

Catalog Number: C10253

Concentration: 1mg/ml

Swiss-Prot No.: O96004

Other Names: Heart- and neural crest derivatives-expressed protein 1, Extraembryonic tissues, heart, autonomic nervous system and neural crest derivatives-expressed protein 1, eHAND, HAND1, EHAND

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

HAND1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: HAND1 antibody detects endogenous levels of total HAND1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 80%, Positives = 80%);
Rat (Identities = 80%, Positives = 80%)

Applications: IF: 1:100~1:500 ELISA: 1:10000

References:

Knoefler M., Gene 224:77-86(1998).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

