



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ACL6A Antibody [#C10255]

Catalog Number: C10255

Concentration: 1mg/ml

Swiss-Prot No.: O96019

Other Names: ACTL6, ACTL6A, MGC5382, actin-like 6A

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ACL6A.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ACL6A antibody detects endogenous levels of total ACL6A protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:


Zhao K., Cell 95:625-636(1998).

Harata M., Biosci. Biotechnol. Biochem. 63:917-923(1999).

Wiemann S., Genome Res. 11:422-435(2001).



Western blot analysis of extracts from Jurkat cells and COLO205 cells, using ACL6A antibody.



Immunofluorescence analysis of HUVEC cells, using ACL6A antibody.