



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MOS Antibody [#C10258]

Catalog Number: C10258

Concentration: 1mg/ml

Swiss-Prot No.: P00540

Other Names: ADRB2, EC 2.7.11.1, Proto-oncogene mos, c-mos, kinase Mos

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MOS.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MOS antibody detects endogenous levels of total MOS protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 IELISA: 1:5000

References:

Watson R., Proc. Natl. Acad. Sci. U.S.A. 79:4078-4082(1982).

The MGC Project Team, Genome Res.

14:2121-2127(2004).

