



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ALDOA Antibody [#C10268]

Catalog Number: C10268

Concentration: 1mg/ml

Swiss-Prot No.: P04075

Other Names: ALDA, ALDOA, ALFA, EC 4.1.2.13,
Fructose-bisphosphate aldolase A, Lung cancer antigen
NY-LU-1, Muscle-type aldolase

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from N-terminal of human
ALDOA.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: ALDOA antibody detects endogenous levels
of total ALDOA protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Sakakibara M., Biochem. Biophys. Res. Commun.
131:413-420(1985).

Izzo P., Eur. J. Biochem. 164:9-13(1987).

Mukai T., Eur. J. Biochem. 195:781-787(1991).

