



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

p53 Antibody [#C10276]

Catalog Number: C10276

Concentration: 1mg/ml

Swiss-Prot No.: P04637

Other Names: Antigen NY-CO-13, Cellular tumor antigen p53, Phosphoprotein p53, TP53, Tumor suppressor p53

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human p53.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: p53 antibody detects endogenous levels of total p53 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Zakut-Houri R., EMBO J. 4:1251-1255(1985).

Lamb P., Mol. Cell. Biol. 6:1379-1385(1986).

Harlow E., Mol. Cell. Biol. 5:1601-1610(1985).



For Research Use Only