



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Histone 1F0 Antibody [#C10293]

Catalog Number: C10293

Concentration: 1mg/ml

Swiss-Prot No.: P07305

Other Names: H1, H1 histone family, member 0, H10, H1FV, histone H1', histone H1.0

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Histone 1F0.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Histone 1F0 antibody detects endogenous levels of total Histone 1F0 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Doenecke D., J. Mol. Biol. 187:461-464(1986).

Collins J.E., Genome Biol.

5:RESEARCH84.1-RESEARCH84.11(2004).

Dunham I., Nature 402:489-495(1999).

