



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

FES Antibody [#C10294]

Catalog Number: C10294

Concentration: 1mg/ml

Swiss-Prot No.: P07332

Other Names: Proto-oncogene tyrosine-protein kinase Fes/Fps, EC 2.7.10.2, C-Fes, FES, FPS

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FES.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FES antibody detects endogenous levels of total FES protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Alcalay M., Oncogene 5:267-275(1990).

Roebroek A.J.M., EMBO J. 4:2897-2903(1985).

Scott A., J. Biomol. NMR 31:357-361(2005).

