



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

HEXB Antibody [#C10295]

Catalog Number: C10295

Concentration: 1mg/ml

Swiss-Prot No.: P07686

Other Names: Beta-hexosaminidase beta chain [Precursor], EC 3.2.1.52, N-acetyl-beta-glucosaminidase, Beta-N-acetylhexosaminidase, Hexosaminidase B, Cervical cancer proto-oncogene 7 protein, HCC-7, Beta-hexosaminidase beta-B chain, Beta-hexosaminidase beta-A chain, HEXB, HCC7

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human

HEXB.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: HEXB antibody detects endogenous levels of total HEXB protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100


ELISA: 1:20000

References:

Korneluk R.G., J. Biol. Chem. 261:8407-8413(1986).


Neote K., Genomics 3:279-286(1988).

Proia R.L., Proc. Natl. Acad. Sci. U.S.A. 85:1883-1887(1988).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using HEXB antibody.



Western blot analysis of extracts from Jurkat cells, using HEXB antibody.