



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

HOXC6 Antibody [#C10318]

Catalog Number: C10318

Concentration: 1mg/ml

Swiss-Prot No.: P09630

Other Names: Homeobox protein Hox-C6, Hox-3C, HHO.C8, CP25, HOXC6, HOX3C

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human HOXC6.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: HOXC6 antibody detects endogenous levels of total HOXC6 protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Simeone A., Proc. Natl. Acad. Sci. U.S.A. 84:4914-4918(1987).

Simeone A., Nucleic Acids Res. 16:5379-5390(1988).

Chariot A., Biochem. J. 319:91-97(1996).



For Research Use Only