



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

FGR Antibody [#C10321]

Catalog Number: C10321

Concentration: 1mg/ml

Swiss-Prot No.: P09769

Other Names: C-FGR, EC 2.7.10.2, FGR protein, P55-FGR, Proto-oncogene tyrosine-protein kinase FGR, SRC2, kinase Fgr

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FGR.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FGR antibody detects endogenous levels of total FGR protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Katamine S., Mol. Cell. Biol. 8:259-266(1988).

Gregory S.G., Nature 441:315-321(2006).

Brickell P.M., Br. J. Cancer 58:704-709(1988).

