



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

MYB-A Antibody [#C10332]

Catalog Number: C10332

Concentration: 1mg/ml

Swiss-Prot No.: P10243

Other Names: AMYB, MYBA, MYBL1, Myb-related protein A

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MYB-A.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MYB-A antibody detects endogenous levels of total MYB-A protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Nomura N., Nucleic Acids Res. 16:11075-11089(1988).

Golay J., Oncogene 9:2469-2479(1994).

Ma X.P., Cancer Res. 54:6512-6516(1994).

