



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

KAP0 Antibody [#C10337]

Catalog Number: C10337

Concentration: 1mg/ml

Swiss-Prot No.: P10644

Other Names: cAMP-dependent protein kinase type I-alpha regulatory subunit, Tissue-specific extinguisher 1, TSE1, PRKAR1A, PKR1, PRKAR1, TSE1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human KAP0.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KAP0 antibody detects endogenous levels of total KAP0 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Sandberg M., Biochem. Biophys. Res. Commun. 149:939-945(1987).

Sandberg M., Biochem. Biophys. Res. Commun. 167:323-330(1990).

Jones K.W., Cell 66:861-872(1991).

