



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

KPSH1 Antibody [#C10346]

Catalog Number: C10346

Concentration: 1mg/ml

Swiss-Prot No.: P11801

Other Names: Serine/threonine-protein kinase H1, EC 2.7.11.1, PSK-H1, PSKH1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human KPSH1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KPSH1 antibody detects endogenous levels of total KPSH1 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Brede G., Genomics 70:82-92(2000).

Hanks S.K., Proc. Natl. Acad. Sci. U.S.A. 84:388-392(1987).

Brede G., Nucleic Acids Res. 30:5301-5309(2002).

