



Web: www.anbio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

HARS Antibody [#C10349]

Catalog Number: C10349

Concentration: 1mg/ml

Swiss-Prot No.: P12081

Other Names: EC 6.1.1.21, HRS, HisRS, histidine transylase, histidine-tRNA ligase, histidyl-tRNA synthetase

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human HARS.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: HARS antibody detects endogenous levels of total HARS protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:

Raben N., Nucleic Acids Res. 20:1075-1081(1992).
Tsui F.W.L., Nucleic Acids Res. 15:3349-3367(1987).
Tsui H.W., Gene 131:201-208(1993).



Western blot analysis of extracts from HeLa cells, using HARS antibody.



Immunofluorescence analysis of HepG2 cells, using HARS antibody.