



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## SKIL Antibody [#C10353]

**Catalog Number:** C10353

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P12757

**Other Names:** Ski-like protein, Ski-related protein,  
Ski-related oncogene, SKIL, SNO

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from C-terminal of human SKIL.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** SKIL antibody detects endogenous levels of total SKIL protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

Nomura N., Nucleic Acids Res. 17:5489-5500(1989).

Pearson-White S.H., Nucleic Acids Res. 25:2930-2937(1997).

Pearson-White S., Nucleic Acids Res. 21:4632-4638(1993).

