



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

UNG Antibody [#C10355]

Catalog Number: C10355

Concentration: 1mg/ml

Swiss-Prot No.: P13051

Other Names: DGU, EC 3.2.2.-, UDG, UNG15, uracil-DNA glycosylase

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human UNG.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: UNG antibody detects endogenous levels of total UNG protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Olsen L.C., EMBO J. 8:3121-3125(1989).

Haug T., FEBS Lett. 353:180-184(1994).

Nilsen H., Nucleic Acids Res. 25:750-755(1997).

