



Web: [www.anbobio.com](http://www.anbobio.com)

[www.dijibio.com](http://www.dijibio.com)

Order: [order@dijibio.com](mailto:order@dijibio.com)

Tel: (86) 519 8805 0026

## TCTP Antibody [#C10362]

**Catalog Number:** C10362

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P13693

**Other Names:** HRF, TCTP, fortilin, histamine-releasing factor, translationally-controlled tumor protein

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human TCTP.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** TCTP antibody detects endogenous levels of total TCTP protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 IHC: 1:50~1:100


ELISA: 1:20000

**References:**

Gross B., Nucleic Acids Res. 17:8367-8367(1989).


Dunham A., Nature 428:522-528(2004).

Hochstrasser D.F., Electrophoresis 13:992-1001(1992).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using TCTP antibody.



Western blot analysis of extracts from HepG2 cells, COLO205 cells and HeLa cells, using TCTP antibody.