



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## HOXA1/B1/D1 Antibody [#C10372]

**Catalog Number:** C10372

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Homeobox protein Hox-B1, Hox-2l, HOXB1, HOX2l

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human HOXA1/B1/D1.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** HOXA1/B1/D1 antibody detects endogenous levels of total HOXA1/B1/D1 protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%);  
Mouse (Identities = 100%, Positives = 100%)


**Applications:** IHC: 1:50~1:100 ELISA: 1:10000

**References:**

Acampora D., Nucleic Acids Res. 17:10385-10402(1989).

Boncinelli E., Genome 31:745-756(1989).

Piper D.E., Cell 96:587-597(1999).



**For Research Use Only**