



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

hnRNP L Antibody [#C10375]

Catalog Number: C10375

Concentration: 1mg/ml

Swiss-Prot No.: P14866

Other Names: HNRPL, HNRPL;, heterogeneous nuclear ribonucleoprotein L, hnRNP-L

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human hnRNP L.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: hnRNP L antibody detects endogenous levels of total hnRNP L protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
IF: 1:100~1:500 ELISA: 1:20000

References:

Pinol-Roma S., J. Cell Biol. 109:2575-2587(1989).


Ito M., Cancer Res. 61:2038-2046(2001).

Rasmussen H.H., Electrophoresis 13:960-969(1992).




Peptide - +

Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using hnRNP L antibody.



Western blot analysis of extracts from HeLa cells, Jurkat cells and HUVEC cells, using hnRNP L antibody.



Peptide - +

Immunofluorescence analysis of HeLa cells, using hnRNP L antibody.