



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Fra-1 Antibody [#C10381]

Catalog Number: C10381

Concentration: 1mg/ml

Swiss-Prot No.: P15407

Other Names: FOSL1, FRA-1, Fos-related antigen 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Fra-1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Fra-1 antibody detects endogenous levels of total Fra-1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Matsui M., *Oncogene* 5:249-255(1990).

The MGC Project Team, *Genome Res.* 14:2121-2127(2004).

Tsuchiya H., *J. Virol.* 67:7001-7007(1993).

