



Web: www.anbio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TCF3 Antibody [#C10385]

Catalog Number: C10385

Concentration: 1mg/ml

Swiss-Prot No.: P15923

Other Names: E12, ITF1, Immunoglobulin enhancer binding factor E12/E47, Immunoglobulin transcription factor-1, Kappa-E2-binding factor, TCF-3, TCF3, TFE2, Transcription factor E2-alpha, Transcription factor ITF-1, Transcription factor-3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human

TCF3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TCF3 antibody detects endogenous levels of total TCF3 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:20000

References:

Kamps M.P., Cell 60:547-555(1990).


Nourse J., Cell 60:535-545(1990).

Grimwood J., Nature 428:529-535(2004).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue using TCF3 antibody.



Peptide - +

Western blot analysis of extracts from 293 cells, using TCF3 antibody.