



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PFKFB1 Antibody [#C10389]

Catalog Number: C10389

Concentration: 1mg/ml

Swiss-Prot No.: P16118

Other Names: 6-phosphofructo-2-kinase/fructose-2,6-biphosphatase 1; 6PF-2-K/Fru- 2, 6-P2ASE liver isozyme; 6-phosphofructo-2-kinase; Fructose-2,6-bisphosphatase; PFKFB1; F6PK, PFRX

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human

PFKFB1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PFKFB1 antibody detects endogenous levels of total PFKFB1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Lange A.J., Nucleic Acids Res. 18:3652-3652(1990).

Ross M.T., Nature 434:325-337(2005).

Algaier J., Biochem. Biophys. Res. Commun. 153:328-333(1988).

