



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## TACD1 Antibody [#C10390]

**Catalog Number:** C10390

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P16422

**Other Names:** Tumor-associated calcium signal transducer 1 [Precursor], Major gastrointestinal tumor-associated protein GA733-2, Epithelial cell surface antigen, Epithelial glycoprotein, EGP,

Adenocarcinoma-associated antigen, KSA, KS 1/4 antigen, Cell surface glycoprotein Trop-1, CD326 antigen, TACSTD1, GA733-2, M1S2, M4S1, MIC18, TROP1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against

synthesized peptide derived from internal of human TACD1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** TACD1 antibody detects endogenous levels of total TACD1 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Strnad J., Cancer Res. 49:314-317(1989).

Simon B., Proc. Natl. Acad. Sci. U.S.A. 87:2755-2759(1990).

Perez M.S., J. Immunol. 142:3662-3667(1989).

