



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

FER Antibody [#C10391]

Catalog Number: C10391

Concentration: 1mg/ml

Swiss-Prot No.: P16591

Other Names: EC 2.7.10.2, FER, FERT2, TYK3, c-FER, kinase Fer, p94-FER

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human FER.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FER antibody detects endogenous levels of total FER protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Hao Q.-L., Mol. Cell. Biol. 9:1587-1593(1989).

Lee S.-T., Oncogene 8:3403-3410(1993).

Krolewski J.J., Oncogene 5:277-282(1990).

