



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## ZNF18 Antibody [#C10394]

**Catalog Number:** C10394

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P17022

**Other Names:** Zinc finger protein 18, Zinc finger protein KOX11, Heart development-specific gene 1 protein, Zinc finger protein with KRAB and SCAN domains 6, ZNF18, HDSG1, KOX11, ZKSCAN6

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

ZNF18.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ZNF18 antibody detects endogenous levels of total ZNF18 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Guo L.L., Yi Chuan 27:523-530(2005).

Ota T., Nat. Genet. 36:40-45(2004).

Thiesen H.-J., New Biol. 2:363-374(1990).

