



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## ZNF23 Antibody [#C10396]

**Catalog Number:** C10396

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P17027

**Other Names:** Zinc finger protein 23, Zinc finger protein 359, Zinc finger protein KOX16, ZNF23, KOX16, ZNF359

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150m M NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human ZNF23.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ZNF23 antibody detects endogenous levels of total ZNF23 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:1000

**References:**

Zhou L., Biochem. Biophys. Res. Commun. 295:862-868(2002).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Thiesen H.-J., New Biol. 2:363-374(1990).

