



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ZNF24 Antibody [#C10397]

Catalog Number: C10397

Concentration: 1mg/ml

Swiss-Prot No.: P17028

Other Names: Zinc finger protein 18, Zinc finger protein KOX11, Heart development-specific gene 1 protein, Zinc finger protein with KRAB and SCAN domains 6, ZNF18, HDSG1, KOX11, ZKSCAN6

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human

ZNF24.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ZNF24 antibody detects endogenous levels of total ZNF24 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Guo L.L., Yi Chuan 27:523-530(2005).

Ota T., Nat. Genet. 36:40-45(2004).

Thiesen H.-J., New Biol. 2:363-374(1990).

