



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

UBF1 Antibody [#C10403]

Catalog Number: C10403

Concentration: 1mg/ml

Swiss-Prot No.: P17480

Other Names: Autoantigen NOR-90, Nucleolar transcription factor 1, TCFUBF, UBF-1, UBF1, UBTF, Upstream binding factor 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human UBF1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: UBF1 antibody detects endogenous levels of total UBF1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:20000


References:

Jantzen H.M., Nature 344:830-836(1990).
Chan E.K.L., J. Exp. Med. 174:1239-1244(1991).
Molina H., Proc. Natl. Acad. Sci. U.S.A. 104:2199-2204(2007).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using UBF1 antibody.



Western blot analysis of extracts from HeLa cells and Jurkat cells, using UBF1 antibody.

For Research Use Only