



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## HOXB9 Antibody [#C10404]

**Catalog Number:** C10404

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P17482

**Other Names:** Homeobox protein Hox-B9, Hox-2E, Hox-2.5, HOXB9, HOX2E

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human HOXB9.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** HOXB9 antibody detects endogenous levels of total HOXB9 protein.

**Reactivity:** Human, Mouse


**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Kosaki K., Teratology 65:50-62(2002).

Giampaolo A., Differentiation 40:191-197(1989).

Boncinelli E., Genome 31:745-756(1989).



**For Research Use Only**