



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## ATF7 Antibody [#C10409]

**Catalog Number:** C10409

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P17544

**Other Names:** ATFA, activating transcription factor 7, cyclic AMP-dependent transcription factor ATF-7

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human ATF7.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ATF7 antibody detects endogenous levels of total ATF7 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Gaire M., Nucleic Acids Res. 18:3467-3473(1990).

Chatton B., Submitted (DEC-1990) to the EMBL/GenBank/DDBJ databases.

Olsen J.V., Cell 127:635-648(2006).

